

## **CHEM525 Experiment 5 Buffer and Gel Preparation**

### **Purpose and Theory:**

Over the next two months, we'll be purifying LDH from *E. coli* cells and then spending several weeks examining its activity under a variety of experimental conditions. The goal of lab this week is to prepare all buffers, gels, and solutions that we will need.

### **Safety Precautions**

- Always wear glove, safety coat and goggles when in the lab.
- If a significant amount of any chemical is spilled, immediately seek the instructor for clean-up protocols.
- Acrylamide is a neurotoxin. Handle it with care. Once polymerized, it is no longer dangerous.

### **Procedure**

1. Make each of the following:
  - a. 1 L of 50 mM Tris Buffer at pH 8.0 with 50 mM NaCl.
  - b. 250 mL of 200 mM Imidazole in the buffer made in part a.
  - c. 250 mL of 1 M NaCl in 50 mM Tris Buffer at pH 8.0.
  - d. 500 mL of 20 mM Tris Buffer at pH 8.0 with 50 mM NaCl.
  - e. 250 mL of 20 mM Tris Buffer at pH 8.5 with 50 mM NaCl.
  - f. 250 mL of Bis-Tris Buffer at pH 6.0 with 50 mM NaCl.
  - g. 250 mL of Bis-Tris Buffer at pH 7.0 with 50 mM NaCl.
2. Make 1 L of Lauria Bertani Broth (25 g/L)
3. Make 1 15% SDS-PAGE Gel according to instructor directions.